



# Cell-IN

enhance cellular uptake

## PROTOCOL

### DESCRIPTION

Cell-IN is a loading reagent based on a polymer formulation manufactured at the Institute of Physical Chemistry Polish Academy of Sciences. Cell-IN ensures effective and cellular uptake of different compounds. Cell-IN is efficient on a wide variety of cell lines.

<b>1. LOADING PROTOCOL FOR 2D CELL CULTURE .....</b>	<b>3</b>
1.1 CELL SEEDING .....	3
1.2 MEDIUM PREPARATION .....	3
1.3 LOADING PROTOCOL .....	3
<b>2. LOADING PROTOCOL FOR CELLS IN SUSPENSION .....</b>	<b>5</b>
2.1 CELL SEEDING .....	5
2.2 MEDIUM PREPARATION .....	5
2.3 LOADING PROTOCOL .....	5
<b>3. TRANSFECTION PROTOCOL.....</b>	<b>7</b>
3.1 CELL SEEDING .....	7
3.2 MEDIUM PREPARATION .....	7
3.3 TRANSFECTION PROTOCOL.....	7
<b>4. TROUBLESHOOTING .....</b>	<b>9</b>
<b>5. PRODUCT INFORMATION .....</b>	<b>10</b>
5.1 ORDERING INFORMATION .....	10
5.2 CONTENT .....	10
5.3 REAGENT USE AND LIMITATIONS.....	10
5.4 QUALITY CONTROL.....	10
5.5 FORMULATION AND STORAGE .....	10
5.6 TRADEMARKS .....	10

# 1. LOADING PROTOCOL FOR 2D CELL CULTURE

## 1.1 Cell seeding

For optimal compound loading conditions, we recommend using cells that are 70 to 80% confluent at the time of loading. Typically, for experiments in 8-well plates, 30 000 cells\* are seeded per well in 300  $\mu$ L of cell growth medium 24h before loading. For other culture formats, refer to Table 1.

\*highly proliferating cells, doubling time  $\sim$ 24h

**Table 1.** Recommended number of cells to seed the day before transfection examples.

Culture vessel*	Number of cells to seed	Surface area per well (cm <sup>2</sup> )	The volume of medium per well to seed the cells (mL)
96-well	8 000 - 10 000	0.3	0.1
24-well	60 000 - 80 000	1.9	0.5
12-well	80 000 - 140 000	3.8	1
8-well	10 000 - 30 000	0.9	0.3
6-well / 35 mm	150 000 - 250 000	9.4	2
60 mm	250 000 - 800 000	21	5

\*for other vessels contact us

## 1.2 Medium preparation

Before compound loading, two media hypertonic and hypotonic are needed to be prepared.

### a. Hypertonic medium

I. Dissolve Cell-IN in serum-free medium in amount according to Table 2

**Table 2.** The recommended amount of medium for dissolving Cell-IN.

Number	Name	The volume of medium (mL)
1	Cell-IN Basic 0.1	0.1
2	Cell-IN Basic 0.5	0.5
3	Cell-IN Basic 1	1
4	Cell-IN Basic 5	5

II. If possible use an ultrasonic bath (36°C) for 30 minutes to enhance the effectiveness of the dissolving process

### b. Hypotonic medium

I. Mix serum-free low glucose medium with sterile water at a rate of 6:4, e.g. 6 mL medium and 4 mL water

## 1.3 Loading protocol

I. Prewarm hypotonic medium and culture medium

II. Prepare working solution by mixing hypertonic medium with your stock solution in a ratio of 9:1, e.g. 18  $\mu$ L hypertonic medium and 2  $\mu$ L stock solution\*

III. Remove medium from cells

IV. Add working solution in amount according to Table 3

\*Please note that to obtain good efficiency the concentration of stock solution should be no less than 100 µM.

**Table 3.** The recommended amount of working solution for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of working solution per well (µL)
96-well	0.3	4
24-well	1.9	24
12-well	3.8	48
8-well	0.9	12
6-well / 35 mm	9.4	125
60 mm	21	280

\*for other vessels, please contact us

V. Incubate for exact 10 minutes in 37°C

VI. Remove working solution

VII. Add hypotonic medium in amount according to Table 4

**Table 4.** The recommended amount of hypotonic medium for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of hypotonic medium per well (µL)
96-well	0.3	90
24-well	1.9	460
12-well	3.8	900
8-well	0.9	280
6-well / 35 mm	9.4	1 800
60 mm	21	4 600

\* for other vessels, please contact us

VIII. Incubate for exact 1.5 minutes in 37°C

CAUTION! Incubation time should not exceed 2 minutes, due to possible cell swelling

IX. Remove hypotonic medium.

X. Add culture medium in amount according to Table 5.

**Table 5.** The recommended amount of culture medium for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of culture medium per well (µL)
96-well	0.3	100
24-well	1.9	500
12-well	3.8	1 000
8-well	0.9	300
6-well / 35 mm	9.4	2 000
60 mm	21	5 000

\* for other vessels, please contact us

XI. Incubate for at least 10 minutes at 37°C.

XII. Change culture medium to medium in which measurements are conducted

## 2. LOADING PROTOCOL FOR CELLS IN SUSPENSION

### 2.1 Cell seeding

For optimal compound loading conditions, we recommend using cells that are 70 to 80% confluent at the time of loading.

### 2.2 Medium preparation

Before compound loading, two media: hypertonic and hypotonic are needed to be prepared.

#### a. Hypertonic medium

I. Dissolve Cell-IN in serum-free medium in amount according to Table 6

**Table 6.** The recommended amount of medium for dissolving Cell-IN.

Number	Name	The volume of medium (mL)
1	Cell-IN Basic 0.1	0.1
2	Cell-IN Basic 0.5	0.5
3	Cell-IN Basic 1	1
4	Cell-IN Basic 5	5

II. If possible use an ultrasonic bath (36°C) for 30 minutes to enhance the effectiveness of the dissolving process

#### b. Hypotonic medium

I. Mix serum-free low glucose medium with sterile water at a rate of 6:4, e.g. 6 mL medium and 4 mL water

### 2.3 Loading protocol

I. Prewarm hypotonic medium and culture medium

II. Prepare working solution by combining hypertonic medium with your stock solution at a rate of 9:1, e.g. 18  $\mu$ L hypertonic medium and 2  $\mu$ L stock solution

III. Transfer 1 mL of medium with cells into a centrifuge tube

IV. Centrifuge medium with cells for 1 minute 2 500 rpm and remove the supernatant

V. Add 20  $\mu$ L of working solution

VI. Incubate for exact 10 minutes in 37°C

VII. Add 1 mL of hypotonic medium

VIII. Incubate for exact 1 minute in 37°C

IX. Centrifuge hypotonic medium for 1 minute 2 500 rpm and remove the supernatant

X. Add 1 mL of culture medium

XI. Incubate for at least 10 minutes in 37°C

XII. Transfer cells into culture vessel and add medium in which measurements are conducted

## 3. TRANSFECTION PROTOCOL

### 3.1 Cell seeding

For optimal compound transfection conditions, we recommend using cells that are 50 to 60% confluent at the time of transfection. Typically, for experiments in 8-well plates, 30 000 cells are seeded per well in 300  $\mu$ L of cell growth medium 24h before transfection and 48h before experiments. For other culture formats, refer to Table 7.

**Table 7.** Recommended number of cells to seed the day before transfection examples.

Culture vessel	Number of cells to seed	Surface area per well (cm <sup>2</sup> )	The volume of medium per well to seed the cells (mL)
96-well	8 000 - 10 000	0.3	0.1
24-well	60 000 - 80 000	1.9	0.5
12-well	80 000 - 140 000	3.8	1
8-well	10 000 - 30 000	0.9	0.3
6-well / 35 mm	150 000 - 250 000	9.4	2
60 mm	250 000 - 800 000	21	5

\* for other vessels, please contact us

### 3.2 Medium preparation

Before compound transfection, two media: hypertonic and hypotonic are needed to be prepared.

#### a. Hypertonic medium

I. Dissolve Cell-IN in serum-free medium in amount according to Table 8

**Table 8.** The recommended amount of medium for dissolving Cell-IN.

Number	Name	The volume of medium (mL)
1	Cell-IN Basic 0.1	0.1
2	Cell-IN Basic 0.5	0.5
3	Cell-IN Basic 1	1
4	Cell-IN Basic 5	5

II. If possible use an ultrasonic bath (36°C) for 30 minutes to enhance the effectiveness of the dissolving process

#### b. Hypotonic medium

I. Mix serum-free low glucose medium with sterile water at a rate of 6:4, e.g. 6 mL medium and 4 mL water

### 3.3 Transfection protocol

I. Prewarm hypotonic medium and culture medium

II. Prepare working solution by mixing hypertonic medium with 125-250 ng/ $\mu$ L\* of desired nucleic acid (NA) in a ratio of 9:1, e.g. 18  $\mu$ L hypertonic medium and 2  $\mu$ L NA

- III. Centrifuge working solution in a vortex centrifuge
- IV. Incubate the working solution at room temperature for 5 min
- V. Remove medium from cells
- VI. Add working solution in amount according to Table 9

\*Please note, that the appropriate concentration of the compound strongly depends on many parameters, such as the size of DNA to be loaded. For the first time, we recommend trying a concentration of 200 ng/ $\mu$ L.

**Table 9.** The recommended amount of working solution for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of working solution per well ( $\mu$ L)
96-well	0.3	4
24-well	1.9	24
12-well	3.8	48
8-well	0.9	12
6-well / 35 mm	9.4	125
60 mm	21	280

\* for other vessels, please contact us

- VII. Incubate for exact 10 minutes in 37°C
- VIII. Remove working solution
- IX. Add hypotonic medium in amount according to Table 10

**Table 10.** The recommended amount of hypotonic medium for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of hypotonic medium per well ( $\mu$ L)
96-well	0.3	90
24-well	1.9	460
12-well	3.8	900
8-well	0.9	280
6-well / 35 mm	9.4	1 800
60 mm	21	4 600

\* for other vessels, please contact us

- X. Incubate for exact 1.5 minutes in 37°C
- CAUTION! Incubation time should not exceed 2 minutes, due to possible cell swelling
- XI. Remove hypotonic medium
- XII. Add culture medium in amount according to Table 12

**Table 12.** The recommended amount of culture medium for the different vessels.

Culture vessel	Surface area per well (cm <sup>2</sup> )	The volume of culture medium per well ( $\mu$ L)
96-well	0.3	100
24-well	1.9	500
12-well	3.8	1 000
8-well	0.9	300
6-well / 35 mm	9.4	2 000
60 mm	21	5 000

\* for other vessels, please contact us

- XIII. Gene expression is detectable after 24h and increases with time



## 4. TROUBLESHOOTING

Observations	Actions
Low efficiency	Use compound with an initial higher concentration
Too high efficiency	Use compound with an initial lower concentration

## 5. PRODUCT INFORMATION

### 5.1 Ordering information

Number	Name	Mass of powder (mg)
1	Cell-IN Basic 0.1	11.6
2	Cell-IN Basic 0.5	58
3	Cell-IN Basic 1	116
4	Cell-IN Basic 5	580

### 5.2 Content

Cell-IN Basic 5 reagent is sufficient to perform up to 50 loadings in 8-well plates or 6 loadings in 6-well plates.

### 5.3 Reagent use and Limitations

For research use only. Not for use in humans.

### 5.4 Quality control

Every batch of Cell-IN reagent is tested by loading TRITC-dextran 75 kDa into HeLa cells.

### 5.5 Formulation and Storage

Cell-IN should be stored at 20°C and is stable for at least one year when stored appropriately. However, dissolved Cell-IN in a serum-free medium can also be stored at 4°C and is stable no longer than one month.

### 5.6 Trademarks

How to cite us: "Cell-IN (IPC PAS, Warsaw, Poland)"